



National Aeronautics and
Space Administration

MICROGRAVITY SCIENCE GLOVEBOX



MSFC Microgravity Science &
Applications Department

Workshop on In-Space Fabrication and Repair

The Microgravity Science Glovebox

July 10, 2003

Charles Baugher

Program Manager

NASA, Marshall Space Flight Center, Science Directorate, Huntsville, AL 35812



National Aeronautics and
Space Administration

MICROGRAVITY SCIENCE GLOVEBOX



MSFC Microgravity Science &
Applications Department



MSG Engineering Unit
Maintained at MSFC



National Aeronautics and
Space Administration

MICROGRAVITY SCIENCE GLOVEBOX

MSG FACT SHEET



MSFC Microgravity Science &
Applications Department

Launched on UF2 in June 2002

Work Volume	906mm Wide x 637mm High 500mm Deep (floor), 385mm Deep (top)
Max. Entrance	406 mm Diameter
Power to Investigation (1000 watts Total)	7 amps @ 28 volts, 4 amps @ 5 volts 2 amps @ +12 volts, 2 amps @ -12 volts 8.3 amps @ 120 volts
Heat Dissipation	800 watts from Cold Plate 200 watts from Air Flow
Video	3 Color, 1 B&W camera, 4 Recorders (Sony, Digital)
Other Resources	Gaseous Nitrogen Vacuum Vent






MICROGRAVITY SCIENCE GLOVEBOX

MSG DATA SYSTEM



National Aeronautics and
Space Administration

MSFC Microgravity Science &
Applications Department





MICROGRAVITY SCIENCE GLOVEBOX

FIRST YEAR OF OPERATIONS



National Aeronautics and
Space Administration

MSFC Microgravity Science &
Applications Department



During its initial year the Microgravity Science Glovebox successfully completed its first materials science experiment (SUBSA) and a series of ESA experiments during the Soyuz 5 Mission. Three experiments are currently on-orbit and in progress.





MICROGRAVITY SCIENCE GLOVEBOX

REPAIR STATUS



National Aeronautics and
Space Administration

MSFC Microgravity Science &
Applications Department

Failure History

November

12 volt Current Spike and 5
volt Converter Failure


February

Repair Unit on Progress

Reoccurrence of 12 volt
Current Spike

Recovery

All Subsystems Disconnected at E-Box and
Individually Reconnected – Fault Disappeared





MICROGRAVITY SCIENCE GLOVEBOX

ON-ORBIT TROUBLE SHOOTING



National Aeronautics and
Space Administration

MSFC Microgravity Science &
Applications Department



**Rack Rotated Forward to Allow
Access to Electronics from Rear**





National Aeronautics and
Space Administration

MICROGRAVITY SCIENCE GLOVEBOX

Example Program



MSFC Microgravity Science &
Applications Department

Science Community

- Instrumentation
- Micro-G Modeling

Workshop

Fabrication and
Repair Testbed

Flight Unit

Engineering Community

- Fabrication Expertise
- Hardware



**Fast and Highly
Visible Program**

2003

2004

2005

2006

2007